


Full-digital swine, ovine ultrasound scanner

accurate reproductive diagnosis on the farm

STANDARD
2
YEAR WARRANTY

Technology Driven
2-Year Warranty



Applications

FarmScan® M30

Visit www.bmv.cc for more information
facebook.com/bmvtechnology | twitter.com/bmvtechnology



Why should you have FarmScan® M30?

Early and independent diagnosis saves time and money. FarmScan® M30 is simple, yet very reliable and resistant, it will be suitable for harsh field conditions.

FarmScan® M30 pregnancy in small livestock can be diagnosed at an early stage, e.g. in pigs on the 19th day, in sheep on the 25th day, especially if the breeder is very experienced.

- 7 hours for continuous operation, meeting a long day scanning needs;
- only 325g, and modern ergonomics designing make it comfort to use;
- advanced full digital core technology to obtain high quality images, ensure accurate pregnancy check at early stage;
- 128 frames permanent storage, without loss of image when powered off;
- high quality 3.8" TFT LCD screen providing a clear, transparent image;
- observation of foetus during pregnancy,
- helps to reduce the cost of running a breeding farm through selection of pregnant female animals and appropriate breeding management.

Superior Imaging Quality



Specification

Imaging mode	Real-time imaging, B,M mode	Probe type	Sector, mechanical
Probe Frequency	2.5/3.5/5.0 MHz	Penetration range	18 cm, 16 cm, 15 cm, 14 cm, 13 cm, 12 cm, 11 cm, 10 cm, 9 cm
Scanning angle	90°	Display	High-resolution TFT-LCD 3.8 inch
Interface	Membrane keyboard	Battery	Li-Ion battery pack, 11.2 V, 2600 mAh
Operating time	7 hours of continuous operation (1 battery)	Charging time	70~90 min
Size	(L)24 cm x (W)12 cm x (H)5 cm	Scanner weight	325 g

4/F, BMV of Yinjin Building, Block 71,
 Baoan Centre District, Shenzhen 518010, China
 TEL/FAX:(0086)755-26564580
 Web: www.bmv.cc Email: export@bmv.cc

BMV
 TECHNOLOGY

*Shaping discovery.
 Improving life.*

Specifications subject change without notice