


Professional ultrasound system

With excellent image quality
for small and large animal pregnancy

STANDARD
2
YEAR WARRANTY

Technology Driven
2-Year Warranty



Applications

FarmScan® M50


Visit www.bmv.cc for more information
facebook.com/bmvtechnology | twitter.com/bmvtechnology



The FarmScan® M50 that makes you a specialist

- Only weight 0.48 kg;
- Excellent imaging quality and ergonomic design for fast and reliable detection;
- 5.8 inch WVGA LCD monitor, wide viewing angle, LCD monitor can turn off when use I-Scan® Video Goggle;
- Support USB ultrasound workstation to transfer images and video to your computer, 450~550 frame cine loop;
- Support leather sheath, designed to be resistant to extremely harsh farm conditions;
- Software & Report for reproductive system, and measurement for distance, area, circumference, volume, angel, heart rate;
- Portability by a shoulder belt giving you total freedom of movements.

Superior Imaging Quality



Benefit on pig farm:

FarmScan® M50 can check whether the pig is pregnant or not as well as 18 days after mating , avoiding the waste of non-pregnant feed. At a rough estimate, this action can directly save about 2000USD feed cost for farm for every 100 head of the sow.

More indirect economic benefits

- Increase the litter size;
- Enhance productivity ;
- Find disease and epidemic situation timely.

Accessories – professional your FarmScan® M50



I-Scan® V1
Video Goggles



mechanical sector probe
(2.5/3.5/5.0 Mhz, adjustable)



Li-ion high-power
-capacity battery



leather protector on farm
(easy-clean by running water)

BMV Technology Co., Ltd

4/F, BMV of Yinjin Building, Block 71,
Baoan Centre District, Shenzhen 518010, China
TEL/FAX:(0086)755-26564580
Web: www.bmv.cc Email: export@bmv.cc



*Shaping discovery.
Improving life.*

Specifications subject change without notice
EN-CA-V1.0

A high-end ultrasound system
particularly good for Pig-Sheep-Goat
fertility and reproduction work.

STANDARD
2
YEAR WARRANTY

Technology Driven
2-Year Warranty

L60RC



PORCINE

- Pregnancy control to save money and feed cost.
- Roughly estimation: US\$2000 feed cost can be saved for every 100 sows.
- Which means our scanner will pay for itself in less than half year, even if you have only 100 sows.

OVINE/CAPRINE

- Pregnancy diagnosis
- Fetal Counting
- Fetal Age
- Fetal viability
- Fetal gender
- Ovarian structures
- Evaluate stage of cycle



Applications

FarmScan® L60 Convex



Visit www.bmv.cc for more information
facebook.com/bmvtechnology | twitter.com/bmvtechnology



Why choose the L60 Convex?

High-End Customers Always Say "NO" to the ultrasound Systems with Mechanical Sector Probes, and Prefer to the Electronic Probes which has:

- Superior Imaging Quality with 96 elements
- More Wider Scan Angle
- Longer Life-time
- No need to do oil-injecting all the time
- No possibility of the air-bubbles which affect image quality in Some Mechanical Probes
- More comfortable for animals without motor-motion



I-Scan® is ideal for BMV ultrasound machines. Every veterinarian deserve to own this. You can see the image in bright sunlight or dark situations when I-Scan® is worn on your head and allows you to see the ultrasound image in front of your eyes.



The sure needs by Pig/Sheep/Goat industry:

Features	Pig-Sheep-Goat-SCAN L60 Convex	Other Sector probe ultrasound systems
Scanning Mode	Electronic Probes	Mechanical Sector Probe
Digital Channel	32	-
Crystal Elements	96	One single element
Scan Angle	160 degrees	90 degrees
Video Goggles	Yes	No

BMV Technology Co., Ltd

4/F, BMV of Yinjin Building, Block 71,
Baoan Centre District, Shenzhen 518010, China
TEL/FAX:(0086)755-26564580
Web: www.bmv.cc Email: export@bmv.cc



*Shaping discovery.
Improving life.*

Specifications subject change without notice
EN-CA-V1.0