



Technology Driven  
2-Year Warranty

# Pad Touch® PT50A VET

Easy and portable pad touch ultrasound designed for professional dog breeders.

## Features

- Superior image quality
- Compact design
- Hi-Res 15" LCD capacitive full touch-screen monitor
- Simplified workflow with focused diagnosis
- Multiple transducers to enable extended diagnosis
- Various animal veterinary measurement packages
- Advanced simple user interface performance
- Digital TGC Preset
- Convenient Handle
- Micro Probe Connector
- External keyboard (option)



Simple



Smart



Immediate result

The Pad Touch® PT50A Vet is a pad touch ultrasound scanner designed specifically for dog breeders to independently confirm pregnancy, estimate the number of puppies before birth(determine number of puppies as early as 23rd day of pregnancy), plan the sale of puppies and observe the pregnancy conveniently and cost effectively. Yes, you can!

## Application

Independently confirm pregnancy, estimate the number of puppies before birth(determine number of puppies as early as 23rd day of pregnancy), estimated date of delivery (EDD).




## Specification

B&W mode: B, M, 2B, B/M, 4B  
 Gray scale: 256  
 Display: 15"LCD monitor  
 Scanning depth: 300mm  
 Cine-loop: 334 frame cine loop  
 memory  
 USB ports: 2

Zoom: Pan zoom for the active image  
 Print: Report print for text&graph&image  
 3 active transducer connector  
 Peripheral ports: Video, 2 USB ports, 1 networking port, 1 power inlet

visit [www.bmv.cc](http://www.bmv.cc) for more information





# True Multi-Specialty Touch Veterinary Ultrasound System

designed for dog breeders

“Pad Touch® PT50A Vet designed for professional dog breeders is an excellent touch pad ultrasound system with a clear image quality. The workflow is very well connected .I was much impressed with the ease and scanning speed offered by it” – Said by owner-operator at “Blessed Bulldogs” and Horse Trainer - Dee Toler.

## Simple Touch, Ease of Use

Intuitive easy-to-use features simplify and increase productivity of examinations, and touch-screen monitor is highly durable for dog breeders. Quick diagnosis report can help independently confirm pregnancy, estimate the number of puppies before birth more easily. This in turn, enables fast and timely treatment for the pregnant dogs.

## Ultimate Accuracy

High quality clinical images strengthen users’ diagnostic confidence thus enabling more accurate diagnoses.

## Fast and Time Saving

Quick boot time and accelerated processes help users attend to more patients in fast-paced clinical environments.

## Compact and Lightweight

The fully functional tablet-based ultrasound system is slim and lightweight, at 3 kg (6.7 lb). Users can easily take the system to patient locations.

## Special Designed small probes

More portable, and can storage the probe inside the system.

※BMV is one of the world's leading researchers, developers and manufacturers of ultrasound across medical and veterinary fields. Founded in 2008, driven by an investment of 12 percent of revenues into R&D, its range of machines now covers vet everything from the lightest and most portable of scanners, to the very latest and most sophisticated in ultrasound technology. BMV also produces vet patient monitor and other veterinary products

[www.bmv.cc](http://www.bmv.cc)



Specification subject to change without prior notice  
©2008-2016 BMV Technology Co., Ltd. All Rights Reserved

**BMV Technology Co., Ltd**

Address: 4/F, BMV of Yinjin Building, Block 71,  
Baoan Centre District, Shenzhen , China  
Tel: + 86-755-26564580; Email: [export@bmv.cc](mailto:export@bmv.cc)

