



BMV Veterinary Systems

PRODUCT CATALOG



www.bmv-vet.com

BMV
TECHNOLOGY

BCL-300 Vet

VET Hematology Analyzer

- ▶ **21 Parameters:** WBC. LYM#, MID#, GR#. LYM%, MID%. GR%, RBC. HGB. HCT, MCV. MCH, MCHC. RDW-CV, RDW-SD. PLT. MPV. PDW, PCT.
- ▶ P-LCC. P-LCR
- ▶ **3 Histograms:** WBC. RBC. PLT
- ▶ **Species:** 12 pre-set species including Canine, Feline, Equine, Bovine, Murine. Swine, Leporid 4 optional species
- ▶ **Sample Volume:** Whole blood. 10 µl
- ▶ **Storage Capacity:** Auto-save, more than 200,000 analysis results available
- ▶ **Throughout:** 60 tests / hour
- ▶ **Data Search Mode:** By Sample ID, Patient ID. Date. Multi-keywords.
- ▶ Loading/extracting data, offline reading
- ▶ **Display:** 10.4 inch colorful touch screen
- ▶ **Operation Mode:** Touch screen, Keyboard/mouse
- ▶ **Quality Control:** L-J, X-R, X-B, etc
- ▶ **Printout:** Thermal printer (defaulted), available for external printers.
- ▶ Optional printing templates and format
- ▶ **Data Communication:** LIS, HL7
- ▶ **Wireless Connection:** Wifi, Bluetooth
- ▶ **Ports:** LAN x 1. RS232 x 1, HDMI x 1, USB Port x 4
- ▶ **Languages:** English, Chinese
- ▶ **Optional Accessories:** Handheld Barcode Scanner

BMV
TECHNOLOGY

bmv-vet.com

BCL-105 Vet

Fluorescent Reader With Bluetooth Printer

- ▶ The evaporator is fully inspected from the factory to ensure that each one is accurate
- ▶ Track-type CO2 absorber
- ▶ Apl valve is closed with one key, two circuits are switched with one key, evaporator is disassembled with one key, intelligent weighing is opened with one key
- ▶ Intelligent exhaust gas recovery, intelligent alarm of weighing machine, avoid overload of filter tank
- ▶ Humanized design, block head in breathing circuit, lockable caster, monitor mounting frame
- ▶ Save space, compact design, integrated installation of monitors, ventilator, E-type oxygen cylinder, oxygen generator and other supporting equipment



The State Scientific and Technological Progress Award
National Standards Maker of Hematology Analyzers and Reagents
National Hematology Analyzers and Reagents Example Manufacturing Park
Excellent Domestic Manufacturer Award by National Health Commission

BCL-120 Vet

Veterinary Point-Of-Care Chemistry Analyzer

- ▶ Three in One: veterinary BP technology, pulse oximetry and temperature.
- ▶ Simple to use: Intuitive, color, touch screen interface.
- ▶ Accurate: BP algorithm was specifically developed for companion animals by experts at a leading vet school.
- ▶ Continuous Monitoring: Monitoring during procedures is simple with user-defined intervals for BP measurement.
- ▶ Portable: Battery-operated device easily moved to different monitoring locations.
- ▶ Rechargeable: Battery recharged when device is plugged into AC power, reducing disposable battery waste.
- ▶ Data transfer to IPAD, PC



Clinical & Laboratory

PRODUCT LINE

BCL-200 Vet

Fluorescence Quantitative Analyzer

- ▶ 3 in 1: Veterinary BP technology + pulse oximetry + temperature. Simple to use: Intuitive, color, touch screen interface.
- ▶ Accurate: BP algorithm was specifically developed for companion animals including veterinary patients.
- ▶ Continuous Monitoring: Monitoring during procedures is simple with user-defined intervals for BP measurement.
- ▶ Portable: Battery-operated device easily moved to different monitoring locations.
- ▶ Rechargeable: Battery recharged when device is plugged into AC power, reducing disposable battery waste.
- ▶ Data transfer to IPAD ,PC



Sample ID: 202002080001
 Pet's Name: Jaguar
 Species:
 Age: 0Y 0M
 Sex:
 Sample Type:

 Item: c-Dimer
 Concentration: >10.00 mg/L
 Test Result:
 Recommendation:
 <=0.55mg/L Normal
 >0.55mg/L

 Test Time: 2020-02-08 22:59:38
 Veterinarian: Dr. MEW

Rapid Test Kit




| | | | |
|---|-----------------------|--|--|
| A | Inflammatory products | Canine C-Reactive Protein (cCRP) Kit, 20 TEST KITS/BOX | Diagnosis and differential diagnosis of acute infectious diseases; Post-operative infection monitoring; Observation of the effect of antibiotics; Disease monitoring and prognosis. |
| | | Feline Serum Amyloid A (fSAA) Rapid Test Kit, 20 TEST KITS/BOX | Serum amyloid A can reach its peak within 4-6 hours after the occurrence of systemic inflammation in the body. After the treatment of the disease, SAA level will significantly drop or return to normal. Continuous detection of SAA can monitor the development of the disease, surgical situation, and efficacy of antibiotics. |
| B | Reproductive products | Canine Progesterone (cProg) Kit, 10 TEST KITS/BOX | Progesterone is the most sensitive clinical indicator for the optimal mating period |

Clinical & Laboratory

PRODUCT LINE

| | | | |
|---------|------------------------------------|---|--|
| C | Viral infectious diseases products | <p>PET INFECTIOUS DISEASE (ANTIGEN DETECTION) RAPID TEST KIT, 20 TEST KITS/BOX . This Kits come in 5 flavors to suit your needs and petcare departments: CPV, C-CV, FPV, FCV, FHV</p> | <p>The kits are used for rapid quantitative detection of IgG antibody against canine Distemper virus, Parvovirus and feline Distemper virus . Vaccination is one of the effective way to prevent infectious diseases in pets. Vaccination can not only reduce the susceptibility to infectious diseases but also improve the specific resistance of cats and dogs to infectious diseases. Successful immunization not only depends on the quality of vaccination, vaccination route and immunization procedures and other external conditions but also depends on the immune response ability of the body. Regular antibody measurement is important to the pet to maintain a healthy state.</p> |
| | | | CPV |
| | | | CCV |
| | | | CCV+CPV |
| | | | CDV |
| | | | FPV |
| | | | FCV |
| | | | FHV |
| | | <p>PET INFECTIOUS DISEASE (IGG ANTIBODY DETECTION) RAPID TEST KIT, 10 TEST KITS/-BOX</p> | <p>The kits are used for rapid quantitative detection of canine Distemper virus, canine Parvovirus, canine Coronavirus and feline Distemper virus. The kits can be used for epidemiological investigation of pet infectious diseases to provide accurate disease information for the owners and doctors, and to take preventive or therapeutic measures.</p> |
| | | | CDV IgG |
| | | | CPV IgG |
| | | | CCV IgG |
| | | | FPV IgG |
| | | | FCV IgG |
| FHV IgG | | | |
| D | More RAPID TEST KIT | 10 TEST KITS/BOX | fPL |
| | | | Canine C/Cys C |
| | | | Feline C/Cys C |
| | | | Canine TT4 |
| | | | Feline TT4 |
| | | | Canine TSH |
| | | | Canine cSDMA |
| | | | Feline fSDMA |



Clinical & Laboratory

PRODUCT LINE

Wechat



Download and intall APP



Cooperative Partners



Academic and Clinical Research



China Agricultural University



South China Agricultural University

BMV Technology Co., Ltd

Office Add: Flat A, 5/F., Unionway Commercial Centre, 283 Queen's Road Central, Hong Kong

Factory Add: BMV of Yinjin Building, Baoan Centre District, Shenzhen 518010, China

Tel: +86 (0755) 85275443 Fax: +86 (0755) 26564580

Contact: Mr. David Luo, CEO

Email: sales@bmv.cc Mobile/WhatsApp: +8613500002887

Website: www.bmv.cc or www.bmv-vet.com